ASPINA - DO M

- EN User manual
- DE Benutzerhandbuch
- FR Manuel de l'utilisateur
- RU Руководство пользователя
- PL Instrukcja obsługi
- SK Návod na použitie
- CS Návod k použití

DENTAL SUCTION EQUIPMENT
DENTALE ABSAUGEINHEIT
L'ASPIRATEUR DENTAIRE
ДЕНТАЛЬНОЕ ОТСАСЫВАЮЩЕЕ УСТРОЙСТВО
ODSYSARKA DENTALNA
DENTÁLNA ODSÁVAČKA
DENTÁLNÍ ODSÁVAČKA

ASPINA - DO M

EKOM spol. s r. o.
Priemyselná 5031/18
SK-921 01 Piešťany
Slovak Republic
tel.: +421 33 7967255
fax: +421 33 7967223

www.ekom.sk
email: ekom@ekom.sk

DATE OF LAST REVISION
DATUM DER LETZTEN ÜBERARBEITUNG
DATE DE MISE A JOUR
ДАТА ПОСЛЕДНЕГО ПЕРЕСМОТРА
DATA WYDANIA
DÁTUM POSLEDNEJ REVÍZIE
DATUM POSLEDNÍ REVIZE

12/2021

NP-DO M-39_12-2021
112000080-000

EN	<u>CONTENTS..... 5</u>
DE	<u>INHALT 31</u>
FR	<u>TABLES DES MATIÈRES 57</u>
RU	<u>ОГЛАВЛЕНИЕ 83</u>
PL	<u>SPIS TREŚCI 110</u>
SK	<u>OBSAH 136</u>
CS	<u>OBSAH 162</u>

CONTENTS

IMPORTANT INFORMATION	6
1. CONFORMITY WITH THE REQUIREMENTS OF THE EUROPEAN UNION.....	6
2. INTENDED USE	6
3. CONTRAINDICATIONS AND SIDE-EFFECTS	6
4. SYMBOLS.....	6
5. WARNINGS	7
6. STORAGE AND TRANSPORT	9
PRODUCT DESCRIPTION	10
7. PRODUCT FUNCTION	10
TECHNICAL DATA.....	13
INSTALLATION.....	17
8. INSTALLATION CONDITIONS	17
9. PNEUMATIC CONNECTION	18
10. ELECTRICAL CONNECTION	18
11. COMMISSIONING	18
12. PNEUMATIC AND ELECTRICAL DIAGRAMS	19
OPERATION.....	21
13. SWITCHING ON THE DENTAL ASPIRATOR.....	21
PRODUCT MAINTENANCE	23
14. MAINTENANCE, CLEANING AND DISINFECTION.....	23
15. LONG-TERM SHUTDOWN.....	28
16. DISPOSAL OF DEVICE	28
TROUBLESHOOTING.....	29
17. REPAIR SERVICE	30
ANNEX	188
23. INSTALLATION RECORD	188

IMPORTANT INFORMATION

1. CONFORMITY WITH THE REQUIREMENTS OF THE EUROPEAN UNION

This product conforms to the requirements of the Medical Device Directive (93/42/EEC), and is safe for the intended use if all safety instructions are followed.

2. INTENDED USE

Dental suction equipment DO M is intended for suction, separation and retention of solid and liquid particles generated during dental procedure.

Risk of injury or equipment damage.

The device must not be used to suction aggressive, flammable and explosive mixtures.

Any other use of the product beyond its intended use is considered an incorrect use. The manufacturer is not liable for any damages or injuries resulting from the incorrect use.

3. CONTRAINDICATIONS AND SIDE-EFFECTS

There are no contraindications or side-effects known.

4. SYMBOLS

The following symbols and marks are used in the User manual, on the device and its packaging:

General warnings

Warning

Danger, electric shock hazard.

Read the operating instructions

Refer to instruction manual

CE – marking

Attention! Hot surface.

Handling mark on package – FRAGILE

Handling mark on package – THIS SIDE UP

Handling mark on package – KEEP DRY

Handling mark on package – TEMPERATURE LIMITATIONS

	Handling mark on package – LIMITED STACKING
	Mark on package – RECYCLABLE MATERIAL
	Connection of protective ground wire
	Fuse
	Danger of biological hazard
	Type B applied part
	Indicates a medical device that has not undergone sterilization
	Sterilizable in a steam sterilizer (autoclave) at specified temperature
	Manufacturer

5. WARNINGS

The product is designed and manufactured to be safe for the user and the surrounding environment when used in the defined manner. Keep the following warnings in mind.

5.1. General warnings

READ THE USER MANUAL CAREFULLY BEFORE USING THE DEVICE AND KEEP IT FOR FUTURE USE!

- The user manual supports the correct installation, operation and maintenance of the product. Careful review of this manual will provide the information necessary for the correct operation of the product for its intended use.
- Original packaging should be retained for the return of the device. Only the original packaging ensures protection of the device during transport. If it is necessary to return the product during the warranty period, the manufacturer is not liable for damages caused by improper packaging.
- This warranty does not cover damages originating from the use of accessories other than those specified or recommended by the manufacturer.
- The manufacturer only guarantees the safety, reliability and function of the device if:
 - installation, new settings, changes, expansion, and repairs are performed by the manufacturer or an organisation authorised by the manufacturer.
 - the product is used pursuant to the user manual.
- The user manual corresponds to the configuration of the product and its compliance with the applicable safety and technical standards at the time of printing. The manufacturer reserves all rights for the protection of its configuration, methods and names.
- This user manual is the original instructions. Translation is performed in accordance with the best available knowledge.

5.2. General safety warnings

The manufacturer designed and manufactured the product to minimise all risks when used correctly for the intended use. The manufacturer considers it its obligation to lay down the following general safety precautions.

- Use and operation of the product must comply with all laws and local regulations valid in the place of use. The operator and user are responsible for following all the appropriate regulations in the interests of performing work safely.
- Only the use of original parts guarantees the safety of operating personnel and the flawless operation of the product itself. Only the accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer should be used.
- The operator must ensure that the device is functioning correctly and safely before every use.
- The user must be familiar with the operation of the device.
- Do not use the product in environments with a risk of explosion.
- The user must inform the supplier immediately if any problem directly related to the operation of the device occurs.

5.3. Safety warnings on protection from electric current

- The device must only be connected to a properly installed, earthed socket.
- Before the product is plugged in, ensure that the mains voltage and frequency stated on the product are in compliance with the values of the mains.
- Prior to putting the device into operation check for any damage to the connected pneumatic lines and electrical wiring. Replace damaged pneumatic lines and electrical wirings immediately.
- Immediately disconnect the product from the mains (remove power cord from the socket) in hazardous situations or in the case of a technical malfunction.
- During all repairs and maintenance, ensure that:
 - the mains plug is removed from the power socket
 - emptied separation vessel and equipment cleaned if possible
- The product shall only be installed by a qualified technician.

6. STORAGE AND TRANSPORT

The manufacturer provides sufficient protection of the appliance when it leaves the factory. Therefore the product is protected against damage during storage and transport.

For the transport, always use an original product packing, if possible.

Transport the dental suction equipment in an upright position

Protect the suction equipment during transport against humidity, impurity or extreme temperatures.

The suction equipment in an original packing can be stored in warm, dry and dustfree rooms.

Keep packaging material if possible. If not, please dispose of the packaging material in an environmentally friendly way. Cardboard can be recycled.

Risk of liquid spillage.

Mobile dental suction equipment may only be transported with emptied separation vessel.

Always empty the separation vessel contents before any transport

Storing or shipping the equipment in any conditions other than those specified below is prohibited.

6.1. Ambient conditions

Products may only be stored and transported in vehicles that are free of any traces of volatile chemicals under the following climactic conditions:

Temperature

–25°C to +55°C

Relative humidity

max. 90% (non-condensing)

PRODUCT DESCRIPTION

Fig. 1: The mobile dental suction equipment ASPINA DO M is built-up on a movable base in a noise damping cabinet. Inside the cabinet there is situated a fan cooled suction pump (9) with the electric distribution and the waste entrapping separation vessel (11). In the lower part – under the box – there is situated a silencer with an output filter (14) and a pre-filter (15) ensuring the air filtration from the suction pump. In the upper – narrowed part of the suction equipment - there is situated a holder for suction tubings (2) equipped with cannulas (1), separation automatics and terminal box with fuses. On the lateral part there is situated a main switch (5), over which are placed the indicators for the network (3) and for the state of filling of the separation vessel (4).

7. PRODUCT FUNCTION

7.1. Function Description


Fig. 1: After switching-on the main switch (5) into the position „I“ the network indicator (3) goes on. Taking the suction tubing (6) off the holder (2) will actuate the suction pump (9) and the vacuum occurs at the suction mouthpiece (1). After repeated putting the suction tubing into the holder, the suction pump turns off. When the separation vessel (11) is filled with waste products, the suction pump turns off and the indicator for the separation vessel (4) filling goes on. Then it is necessary to return the suction tubing back into the holder and to empty the separation vessel. During a longer work, mainly with the desalivating mouthpiece, the box temperature may increase; the cooling fan is then automatically actuated. The fan turns off automatically, when the box temperature drops.

7.2. Detailed Description of the Suction Part Function

Fig. 1: The vacuum air flows together with sucked waste products from the oral cavity through the tubing system from the suction mouthpiece (1) at first through the inlet sieve (7), wherein the solid impurities are entrapped. The vacuum air, together with the sucked waste products, free of solid impurities greater than 2 mm, then flows into the separation vessel (11) wherein it is separated from the vacuum air and entrapped into the separation vessel.

Suction air flows through the filter (18) and into the suction unit (9), from which it is pushed through a noise muffler. Therein the output air passes through the output pre-filter (15) and the bacteriologic output filter (14). After passing the filters, the air free of impurities is blown off into a free space under the dental suction equipment.

Fig. 1: Mobile dental suction equipment ASPINA DO M

Descriptions for fig. 1

- | | |
|---------------------------------------|--------------------------------|
| 1 Cannulas | 10 Separator cover |
| 2 Holder | 11 Separation vessel |
| 3 Network indicator | 12 Bowl |
| 4 Indicator for the separation vessel | 13 Filter holder |
| 5 Main switch | 14 Bacteriologic output filter |
| 6 Suction tubing | 15 Output pre-filter |
| 7 Inlet sieve | 16 Fixing lugs |
| 8 Holder | 17 Antifoaming tables |
| 9 Suction pump | 18 Filter |

TECHNICAL DATA

Products are designed to operate in dry, ventilated and indoor dust-free rooms with the following climatic conditions:

Temperature

+5°C to +40°C

Relative humidity

max. 70%

		DO M
Rated voltage / frequency	V / Hz	230 / 50; 230 / 60*; 110 / 60*
Rated current	A	2.6
Maximum flowrate :	l/min	1100
Maximum vacuum	kPa	12
Sound level	dB(A)	≤48
Mode of operation of the suction pump		permanent S 1
Dimensions of the suction equipment w x d x h	mm	354x575x845
Weight of the suction equipment	kg	36
Classification pursuant to STN EN 60 601-1 (IEC 60601-1)		class I
Applied parts type pursuant to STN EN 60 601-1 (IEC 60601-1)		B
Classification pursuant to MDD 93/42/EEC, 2007/47/EC		Ila

(*) - The rated voltage and frequency can be changed by special order

Electromagnetic compatibility declaration

Medical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided below.

Guidance and manufacturer's declaration - electromagnetic emissions		
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal functions. Therefore, the RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	The equipment is not likely to cause any flicker, as the current flow is approx. constant after the start up.	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration - electromagnetic immunity

Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances

The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV 100 kHz repetition frequency Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV L-N ±2 kV L-PE; N-PE Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption, and voltage variations on power supply input lines IEC 60601-4-11	U _T =0%, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) U _T =0%, 1 cycle U _T =70% 25/30 cycles (at 0°) U _T =0%, 250/300 cycles	U _T >=95%, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) U _T >=95%, 1 cycle U _T =70% (30% dipp U _T), 25(50Hz)/30(60Hz) cycles (at 0°) U _T >=95%, 250(50Hz)/300(60Hz) cycle	Mains power quality should be that of a typical commercial or hospital environment. The device stops and restarts automatically at each dip. In this case does not occur unacceptable pressure drop.

NOTE: U_T is the A.C. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity

Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances

The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$, 80 MHz to 800 MHz $d=2,3\sqrt{P}$, 800 MHz to 2,7 GHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Proximity fields from RF wireless communications equipment IEC 61000-4-3	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.

^b Over the frequency range 150 kHz to 80 kHz, field strengths should be less than 3 V/m.

INSTALLATION

Risk of incorrect installation.

Only a qualified professional can install the equipment and place it into operation for the first time. This professional is obliged to train operating staff as to the use and maintenance of the device. Installation and training of all operators shall be confirmed by the installer's signature on the certificate of installation.

Risk of incorrect installation.

ANY MODIFICATION OF THIS EQUIPMENT IS FORBIDDEN!

Risk of damage to the device.

This equipment cannot be used nearby other instruments. If this equipment is used nearby other instruments, the equipment must be observed in order to verify normal operations in the configuration it will be used.

Instruments may be affected electro-magnetically!

You may notice a "new product" odour when you first place the product into service (for a short period of time). This odour is temporary and does not impede the normal use of the product. Ensure the space is properly ventilated after installation.

8. INSTALLATION CONDITIONS

- The equipment may only be installed and operated in dry, well ventilated and dust-free rooms.
- The mobile dental suction equipment should be installed with regard to an easy access for the operation and maintenance and good accessible rating plate.
- The equipment should stand on a plane, sufficiently stable base (be aware of the suction equipment weight, see Technical data).
- The equipment may not be operated outdoors or in wet or damp environments.

Risk of explosion.

Do not use the equipment in the presence of explosive gases, dust or combustible liquids.

- To guarantee a trouble-free operation of the suction equipment, the room temperature should never drop under +5°C and exceed +40°C. The ideal ambient temperatures are between +10°C up to +25°C.
- Most of the electric energy used by the suction pump (9) is converted to heat and given off to the ambient. During a longer lasting work, above all with the desalivating mouthpiece, the box temperature increases over 40°C causing an automatic actuation of the cooling fan. After the room is cooled down under approximately 32°C, the fan turns off again.

9. PNEUMATIC CONNECTION

Risk of damage to pneumatic components.

Suction tubings must not be broken.

10. ELECTRICAL CONNECTION

- The dental suction equipment model 230V is delivered with the safety plug.
- Keep the socket easily accessible to ensure that in an emergency, the device can be safely disconnected from the mains.

Risk of electric shock.

It is necessary to follow all local electro technical regulations. The mains voltage and frequency must comply with the data stated on the device label.

Risk of fire and electric shock.

Power cord must not be broken.

11. COMMISSIONING

After unpacking, put the suction equipment on the floor, open the door and check whether the separator's cover (10) with sensing electrodes is attached to the separation vessel (11). In case of need, fix the cover to the vessel and then follow the instructions in the Chapter No. 13. Close the door and connect the suction equipment power cord to the socket. Change the main switch (5) over into the position "I", the indicating lamp for the network (3) goes on and the suction equipment is ready for the operation.

Risk of use of incorrect cannulas.

Suction cannulas must comply with local legislative and technical requirements for dental suction cannulas and their connection dimensions must be compatible with suction cannulas connectors of 11mm and 16mm diameter.

Cannulas are not delivered in a sterile state!

Suction cannulas must undergo complete preparation cycles before first use and after each patient or use following instructions in Chapter 14.8


Dental suction equipment is not equipped with a backup power supply!

12. PNEUMATIC AND ELECTRICAL DIAGRAMS

12.1. Pneumatic diagram

Aspina DO M

Description to pneumatic diagram:

- | | | | |
|---|-------------------|---|-------------------|
| 1 | suction cannulas | 5 | suction pump |
| 2 | vacuum valve | 6 | output pre-filter |
| 3 | inlet sieve | 7 | output filter |
| 4 | separation vessel | 8 | noise muffler |


12.2. Electrical diagram

Aspina DO M

1/N/PE ~ 230V/50Hz, 230V/60Hz

ELEKTRIC OBJECT OF 1st. CAT.

Type B

Description to electrical diagrams

M	Suction pump 230 V/50 Hz
EV	Fan 230 V/50-60 Hz
HL1,2	Indicator glow lamps
PS	Control electronics
ST	Temperature switch
Q	Power-supply switch
A	Cover of bracket
B	Complete bracket
C	Complete cover
D	Cabinet cover
E	Main board
F	Doors of cabinet

OPERATION

The equipment may only be operated by the trained staff!

Risk of electric shock.

Disconnect from the mains in case of a danger (unplug) !

Risk of burns or fire.

Hot surfaces of the suction pump. Do not touch - Danger of burn.

13. SWITCHING ON AND OFF THE DENTAL SUCTION EQUIPMENT

Switch on the dental suction equipment by turning the main switch (5) in the position „I“.

The suction equipment is put into operation automatically by taking the suction tubing (6) with a mouthpiece (1) off the holder (2). It works uninterruptedly until both tubings are put on the holders again or until the separation vessel (11) is full (during an uniform operation it is filled within about 6 – 10 hours). Suction tips holder is equipped with vacuum regulation of suction equipment depends on dentist needs. The state of filling of the separation vessel is evaluated by the separation automatics which interrupts the suction by the suction pump (9) and filling of the separation vessel is indicated automatically with the indicator lamp (4) going on. At this point it is necessary to put the tubings with suction mouthpieces into the holder and switch the power-supply switch off. Afterwards open the door by pulling the grips situated on the side walls, uncouple rubber fixing lugs from the separator cover (10) and pull the separator's vessel out. Hang the separator's cover into the holder on the left side (8).

Empty the separation vessel content into the sewer, rinse the vessel with water, connect with the separator's cover in a reversal order. Set the separator into a bowl (12) in the suction equipment cabinet (the bowl is detachable). Check the connection of the cover with the vessel, seating of the separation vessel in the bowl and close the door.

Turn off the dental suction equipment by turning the main switch (5) in the position „O“.

13.1. Use of antifoaming tablets

Under specific exhaust conditions, increased exhaust condensate foaming may occur due to the switching off of the exhaust unit when the separation tank is not full. Antifoaming tables (17) must be inserted in order to prevent this from occurring (see the basic accessories); add these tables to the inlet sieve (7). The tablets will then gradually dissolve in the input screen and will significantly decrease exhaust condensate foaming; the tablets also function as disinfection agents.

13.2. Filter use in the separation vessel

The filter in the separation vessel is an additional layer of protection for the suction equipment to prevent moisture from the foam from entering the suction equipment. Foam is generated in the separation vessel under specific suction conditions. Foam-blocking tablets are used to eliminate such foam, see Chapter 13.1

Schedule for adding tablets: Insert 1 or 2 antifoaming tablets into the inlet sieve (7) at the end of work and after cleaning and disinfecting the equipment.

Biohazard warning.

Suction cannulas must undergo complete preparation cycles before first use and after each patient or use following instructions in Chapter 14.8

Used suction cannulas intended for repeated use must be cleaned and steam sterilized after every patient (see Chapter 14.8).

The cannulas are to be exchanged after maximum of 100 cycles of steam sterilizing.

Cannulas must be visually inspected for any signs of wear or tear before they are used. Never use the cannulas if it is damaged: discard it immediately.

Suction cannulas are considered applied parts type B (according to IEC 60601-1).

Risk of overheating of the equipment .

It's forbidden to cover ventilation openings placed on upper sides of equipment.

Danger of injury.

Before switching the equipment on it is necessary to check if door of cabinet is closed.

PRODUCT MAINTENANCE**14. MAINTENANCE, CLEANING AND DISINFECTION**

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used..

The device has been designed and manufactured to keep its maintenance to a minimum. For correct and reliable operation perform the following operations.

Risk of servicing by persons without required qualification.

Repair works exceeding the scope of a general maintenance may only be carried out by a qualified specialist or by the manufacturer's customer service.

Use only spare parts and accessories permitted by the manufacturer.

Risk of electric shock.

Before making any maintenance, repair or cleaning, it is unavoidable to disconnect the equipment from the network (unplug).

Burn hazard.

Allow the equipment to cool down before maintenance and service!

Reconnect the grounding lead if removed during service work to its previous location once the work is complete.

Biohazard warning.

During the contact with the internal parts of aspirator including preparation phase for suction cannulas, always wear EC marked personal protective equipment (PPE).

Empty the separation vessel and clean the equipment before any maintenance or service.

14.1. Maintenance intervals

Time interval	Once a day	once every 15 days	as needed	max. 100 steam cleaning cycles	once every 3 months	Once a year	once every 18 months	Once every 2 years	Chap. replacement	Set of replacement parts	Perform ed by
Use of antifoaming tablets			x						13.1	-	operator
Inlet Sieve Cleaning	x								14.2	-	operator
Lubricating seals and sliding clamps		x							14.3	-	operator
Disinfection of Tubings and Separation Vessel	x								14.4	-	operator
Cleaning and disinfection of the exterior surfaces of the product			x						14.5	-	operator
Replacement of output pre-filter					x				14.6	-	operator
Replacement of output filter						x			14.7	-	operator
Exchange of suction cannula				x					14.8	-	operator
Replacing the filter in the separation vessel							x		14.9	-	operator
Perform "Repeated Test" according to EN 62353								x	14	-	qualified technician

Following works must be carried out from the hygiene and the correct suction equipment function points of view.

14.2. Inlet Sieve Cleaning

Solid particles (sucked together with a liquid component and the air) are entrapped in the inlet sieve (7) during the equipment operation which therefore must be regularly cleaned, by section 14.1 (though, always at the end of each shift).

The inlet sieve cleaning is possible only after the moisture inside suction tubings (6) and a sieve is reduced by means of the several seconds lasting air suction through tubings taken off the holder (2). Then turn the power-supply switch (5) off in the position „O“ and slightly lift the cover of the inlet sieve. Afterwards pull the sieve out of its place gripping the holder, remove solid particles and clean the sieve. Then put the sieve on the previous place into the holder and carry out the assembly in the contrary manner to the dismantling.

Danger of environmental pollution.

If the suction of amalgam particles is expected, it is necessary to empty the inlet sieve content into a closable vessel and forward it to the collecting centres according to local regulations.

14.3. Lubricating seals and sliding closures

O-rings (sealing rings) and sliding terminal clamps (Fig. 1) must be lubricated using a suitable silicone oil for dental applications (e.g. Lubri-Jet spray) every 15 working days.

14.4. Disinfection of Tubings and Separation Vessel

For the separation vessel (11) disinfection it is necessary to empty the same, rinse with warm water and clean mechanically with warm water and a disinfectant.

The disinfection of internal suction tubings and of the separation vessel should be performed by section 14.1, though always at the end of a shift by suction of warm water containing a disinfectant having a limited foaming power through both suction tubings (6) separately. The used disinfectant must be approved in compliance with a valid national legislation on usage for disinfected surface and a type of material.

For the disinfectant use, follow the manufacturer's instructions.

14.5. Cleaning and disinfection of the exterior surfaces of the product

Use neutral agents for cleaning and disinfection of the external surfaces of the product.

The use of aggressive cleaning agents and disinfectants comprising alcohol and chlorides may lead to the damage of the surface and the discoloration of the product.

To find out whether the suction equipment works correctly, the following works should be carried out in the defined maintenance intervals:

14.6. Replacement of the Output Pre-filter

It is necessary to replace the output pre-filter (15) by section 14.1. Before making any replacements turn the suction equipment power-supply switch (5) off. The pre-filter together with the filter (14) are situated on the left lower part of the suction equipment's base.

For its dismantling, catch the filter holder (13) with one hand and with the other hand slew the fixing lugs (16) by 90° by simultaneous pulling down. Remove the filter holder together with the output

filter. Then release the output pre-filter from the cavity wherein also the output filter was situated (check also the state of pollution of the output filter).

For the reassembly, set the output filter into the filter holder, put the output pre-filter onto the filter (direct the output pre-filter with the gummed cloth towards the filter) and insert all the assembly into the cavity space from below. Slew the fixing lugs back by 90° so that they engage into the recessed parts of the filter holder.

14.7. Replacement of Output Filter

During the regular use of the device, the output filter (14) must be replaced according to the interval in Chapter 14.1. The procedure of disassembly and assembly of the filter is the same as for the replacement of output prefilter (15).

Risk of biological contamination.

Neglecting of replacement interval of the output filter (HEPA) may result in contamination of the ambient air.

14.8. Cannulas handling and preparation operation

The preparation instructions have been drafted by supplier of cannulas in accordance with EN ISO 17664 standard requirements.

Danger of biological hazard.

During the preparation phase for suction cannulas, always wear EC marked personal protective equipment (PPE).

TYPE MANUFACTURER OF CANNULAS DELIVERED FOR THE PRODUCT	Type: Monoart aspirator tips EM21 22910103, Monoart aspirator tips EM21 EVO 22931001 Euronda, Via Chizzalunga 1, 36066 Sandrigo, VI, Italia Tel.(+39) 0444 656185 www.euronda.com
WARNINGS	• Store brand new or unused cannulas in a dry, clean and protected environment • Only use EC marked detergents and disinfectants, suitable for use with polypropylene (PP), following the indications concerning concentration, temperature and time of action provided in the manufacturer's instructions. • Collect and clean the cannulas immediately after use. Dried or sticky residue can make cleaning operations difficult or ineffective. • Never use metal brushes to clean such products. • Damaged cannulas must not be reused.
TRANSPORTATION AFTER USE	Place the product in a specific container with disinfectant solution according to the concentrations and soak times indicated by the manufacturer.
PREPARATION FOR CLEANING	If necessary, remove visible dirt from the product and then rinse under running water.
MANUAL CLEANING	• Use a nylon brush and detergent/disinfectant, suitable for use with polypropylene (PP), to clean the product until all dirt has been removed. Use a bottle brush to clean the insides. • Rinse thoroughly under running water (of good microbiological quality) until all detergent/disinfectant residue has been removed. • Dry well using compressed air.
ULTRASONIC CLEANING	• Place the product in a perforated basket and place it inside the ultrasonic bath filled with detergent or a disinfectant solution suitable for polypropylene (PP).

	• Set the wash cycle following the manufacturer's instructions; it is not recommended to set the temperature above 45°C. • Rinse thoroughly under running water (of good microbiological quality) until all detergent residue has been removed. • Check there is no dirt residue on the product; repeat the cleaning cycle whenever necessary. • Dry well using compressed air. Note: ultrasonic cleaning can be used in conjunction with manual and automatic cleaning and disinfection.
THERMAL DISINFECTOR CLEANING	Use cleaning and disinfectant equipment that is EN ISO 15883 standard-compliant. • Position the product using the supplied holders so that the water flows easily and rinses the inner surfaces. • Follow the manufacturer's instructions concerning selection of cleaning/disinfection cycles and detergents to be used. • On completing the cycle, check there is no dirt residue on the product; repeat the cleaning cycle whenever necessary.
OPERATING CHECKS AND INSPECTION	Visually inspect the condition and cleanliness level of the product. Discard all damaged or worn products.
PACKAGING	Select EC marked packaging that meets the requirements of Council Directive 93/42/CEE and addenda, suitable for the product and the sterilization procedure, leaving enough space to prevent tension when sealed.
STEAM STERILIZATION	Fractionated pre-vacuum steam sterilization using EN 13060 or EN 285 standard compliant equipment and validated sterilization procedure in accordance with EN ISO 17665-1 standard requirements. • Use type B cycles. • Sterilization temperature: 134°C. • Stay time: 5 minutes. • Dry time: 10 minutes.
MAXIMUM DURATION	100 cycles
STORAGE	Storage times shall depend on the type of packaging used and the storage environment. Packaging must be protected against dust, humidity and all contamination risks.

14.9. Replacing the filter in the separation vessel


Fig. 2: Replacing the filter in the separation vessel

15. LONG-TERM SHUTDOWN

When the dental suction equipment is not used for a longer time, it is necessary to carry out the cleaning and disinfection of all parts as described in the Articles No. 14.2, 14.4 and then turn the power-supply switch (5) on in the position "I", take the suction tubings (6) off the holder (2) and let the air intake flow through (for approximately 15-20 min) so that the suction system is fully dried out. Afterwards put the tubings into the holder and turn off the power-supply switch in the position "0", disconnect the equipment from the network (unplug) and remove the suction mouthpieces (1) from suction tubings.

16. DISPOSAL OF DEVICE

- Disconnect the appliance from mains.
- To follow the rules of personal hygiene for works with contaminated material
- Clean device according, section 14.
- To separate, label, packing and providing for decontamination of contaminated parts by course of national regulations.
- Dispose of the product following all locally applicable regulations.

Danger of biological hazard and environmental pollution.

Internal parts of the suction equipment may be contaminated with biological material. Hand over the device to a specialised organisation for decontamination before separating and disposing of the unit.

TROUBLESHOOTING

Risk of electric shock.

Before interfering with the equipment, first disconnect it from the mains (remove the power socket).

Troubleshooting may only be performed by a qualified service technician.

Danger of biological hazard and environmental pollution.

Follow the rules of personal hygiene for manipulation with contaminated material. Separate, label, pack and ensure decontamination of contaminated parts according to local regulations.

In case of repair of parts of device which might be contaminated please follow bellow mentioned instruction:

- Disconnect the appliance from mains
- To follow the rules of personal hygiene for works with contaminated material
- Clean device according, section 14.
- To separate, label, packing and providing for decontamination of contaminated parts by course of national regulations
- To perform the repair of damaged parts

Failure	Problem and possible reasons	Way of removal
Suction equipment does not work	Mains voltage absent	Check mains voltage in the socket, activate a protective element in the electric distribution (fuse, circuit breaker)
	Interrupted supply of el. energy	
	Power-supply switch off	Check the power-supply switch position („I“), switch on the power-supply switch (network indicator goes on)
	Interrupted power cord	Replace the defective part
	Loose clamp of terminal board	Tighten the clamps
	Interrupted motor wiring, defective thermal protection	Replace suction pump
	Defective mains fuse of the equipment	Check the state of equipment mains fuse, replace defective fuse (fuses 2 x T 6.3 A are situated in the upper, narrowed part of the suction equipment)
	Full separation vessel	Empty the separation vessel content
	High current consumption	Defective starting capacitor, replace the capacitor
	Suction equipment is/was awash by liquid	Suction equipment pass on repair to service – keep all tools and norms for work with

		contamination material
Suction equipment is switching without taking suction tubings off	Defective control automatics	Replace the defective part
	Defective microswitch in the suction tubings holder	Replace the defective part
Suction equipment is noisy	Cabinet door improper closed	Close the box door
	Defective bearings of the suction pump	Replace the defective bearing for an identical type
	Incorrectly mounted filter holder	Mount the filter holder into a correct position
	Incorrectly mounted cover of the separation vessel	Mount the separation vessel cover into a correct position
	Suction equipment is/was overflowen by liquid.	Suction equipment pass on repair to service – keep all rools and norms for work with contamination material
Suction equipment efficiency reduced	Heavy polluted output pre-filter	Replace the pre-filter (see Maintenance intervals art.No.14.1)
	Heavy polluted output filter	Replace the filter (see Maintenance intervals art.No. 14.1)
	Leaks in the suction system	Check connections, seal untight connections
	Defective control valve in the suction tubings holder	Replace the defective part
	Heavy polluted inlet sieve	Clean the inlet sieve

17. REPAIR SERVICE

Warranty and post-warranty repairs must be done by the manufacturer, its authorized representative, or service personnel approved by the supplier.

Attention.

The manufacturer reserves the right to make changes to the equipment without notice. Any changes made will not affect the functional properties of the equipment..

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1. Product: (model) DO M		2. Serial number:	
3.1. User's name:			
3.2. Address of installation:			
4. Equipment connected to the compressor:			
5. Installation / Commissioning:		6. Contents of operator training:	
Product completeness check **	Y	Description of the product and functions**	Y
	N		N
Documentation completeness check **	Y	Product operation: turning on/off, controls, control procedures, data on the display panel, alarms, operation in alarm conditions**	Y
	N		N
Installation/connection to equipment **	Y	Product maintenance: maintenance intervals, maintenance procedure, service intervals, operating activities**	Y
	N		N
Functional test **	Y	Safety measures, warnings – their meaning and compliance **	Y
	N		N
Notes:			
7. Operator instructed on safety measures, operations and maintenance:			
Name:		Signature:	
Name:		Signature:	
Name:		Signature:	
8. Installation and instruction performed by: First name/Last name		Signature:	
Company:		Address:	
Phone:			
Email:		Date:	
9. Distributor:			
Company:		Address:	
Contact person:			
Phone:		Email:	

** mark with an "X" in points 5 and 6 (Y - yes /N - no). Enter any observations from points 5 and 6 into the "Notes" section

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