DK50 2X4VR/110

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IMPORTANT INFORMATION

1. CE MARKING

Products labeled with the CE mark of compliance meet the safety guidelines (93/42/EEC) of the European Union.

2. WARNINGS

2.1. General warnings

- This Installation, Operation and Maintenance Manual is a part of the appliance and must be kept with the compressor. Careful review of this manual will provide the information necessary for correct operation of the appliance.
- The safety of operating personnel and trouble-free operation of the appliance are guaranteed only if original parts are used. Only accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer can be used.
- If any other accessories or consumable materials are used, the manufacturer cannot be held responsible for the safe operation of the appliance. This guarantee does not cover damages originating from the use of accessories or consumable material other than those specified or suggested by the manufacturer.
- The manufacturer guarantees the safety, reliability and function of the appliance only if:
 - Installation, new settings, amendments, extensions and repairs are performed by the manufacturer or its representative, or a service provider authorized by the manufacturer
 - The appliance is used in accordance with this Installation, Operation and Maintenance Manual
- The manufacturer reserves all rights for the protection of its wiring diagrams, methods and names.
- Translation of Manual for Installation, Operation and Maintenance is carried out in accordance with the best knowledge. In the case of ambiguities, the Slovak version of the text prevails.

2.2. General safety warnings

The manufacturer developed and designed the equipment in such a way so that any risks were excluded if it is used according to intention. The manufacturer considers it to be its obligation to describe the following safety measures in order to exclude residual damages.

- Operation of the appliance must be in compliance with all local codes and regulations.
- Original packaging should be kept for the return of the appliance. Only the original packaging ensures protection of the appliance during transport. If it is necessary to return the appliance during the guarantee period, the manufacturer is not liable for damages caused by improper packaging.
- Each time the appliance is used, the operator must make sure that it is functioning correctly and safely.
- The user must fully understand the operation of the appliance.
- The product is not intended for operation in areas with a risk of explosion.
- If any problem occurs during use of the appliance, the user must inform his supplier immediately.

2.3. Electrical system safety warnings

- The appliance must be connected to earth (grounded).
- Before the appliance is plugged in, make sure that the mains voltage and mains frequency stated on the appliance are the same as the power mains.
- Prior to putting into operation it is necessary to check for possible damage of the equipment and connected air and electric distributions. Damaged pneumatic and electric lines must be immediately replaced.
- Immediately disconnect the appliance from the mains (pull out mains plug) if a technical failure occurs.
- During repairs and maintenance, ensure that:
 - The mains plug is pulled out from the socket
 - Pressure pipes are vented and pressure is released from the air tank.
- The appliance must be installed by an approved, qualified technician.

3. ALERT NOTICES AND SYMBOLS

In the Installation, Operation and Maintenance Manual and on packaging and product, the following labels or symbols are used for important information:

	Information, instructions and cautions for the prevention of damage to health or materials
	Caution! Dangerous electrical voltage
	Read the user manual!
	CE mark of compliance
	Compressor is remote-controlled and may start without warning
	Caution! Hot surface
	Earth (ground) connection
	Terminal for ground connection
	Fuse
	Alternating current
	Handling mark on package – FRAGILE
	Handling mark on package – THIS SIDE UP
	Handling mark on package – KEEP DRY
	Handling mark on package – TEMPERATURE LIMITATIONS
	Handling mark on package – LIMITED STACKING
	Mark on package – RECYCLABLE MATERIAL

4. STORAGE AND TRANSPORT

The compressor is shipped in cardboard that protects the appliance from damage during transport.

Caution! For transport, always use the original packaging and secure the compressor in the upright position.

Protect the compressor from humidity and extreme temperatures during transport and storage. A compressor in its original packaging can be stored in a warm, dry and dust-free area. Do not store near any chemical substances.

Keep packaging material if possible. If not, please dispose of the packaging material in an environmentally friendly way and recycle if possible.

Caution! Before moving or transporting the compressor, release all the air pressure from the tank and hoses and drain the condensed water.

5. TECHNICAL DATA

			DK50 2x4VR/110	DK50 2x4VR/110S
Rated voltage / frequency	V / Hz		3x400/50	3x400/50
Compressor output at pressure of 6 bar	Lit.min ⁻¹		520	520
Compressor output at 6 bar with dryer	MD	Lit.min ⁻¹	450	450
	NDL		400	400
Maximum compressor current	A		9.4	10.2
Maximum compressor current with dryer	A		10.1	10.6
Motor output	kW		2x2.2	2x2.2
Air tank capacity	Lit.		110	110
Working pressure	bar		6.0 – 8.0	6.0 – 8.0
Maximum operating pressure of safety valve	bar		9.0	9.0
Noise	L _{pfA} [dB]		77	58
Mode of operation			Continuous S1	Continuous S1
Mode of operation with dryer			Continuous S1	Continuous S1
Compressor dimensions	W x L x H	mm	1090x500x830	1260x795x1025
Dimensions with dryer W x L x H	MD	mm	1210x670x830	1550x795x1240
	NDL		1210x770x1240	1550x795x1240
Weight	kg		136	218
Weight with dryer	MD	kg	158	240
	NDL		172	252
Air dryer performance				
Atmospheric dew point	MD	°C	-20	
Atmospheric dew point				
Pressure dew point	NDL		- 35	
Configuration pursuant to EN 60 601-1			Type B, class I.	

Climatic conditions during storage and transport
 Temperature : -25°C to +55°C, 24 h to +70°C
 Relative air humidity : 10% to 90 % (no condensation)

Climatic operation conditions
 Temperature : +5°C to +40°C
 Relative air humidity : 70%

5.1. FAD efficiency correction for differences in elevation

FAD correction table

Elevation [mamsl]	0 - 1500	1501 - 2500	2501 - 3500	3501 - 4500
FAD [l/min]	FAD x 1	FAD x 0.8	FAD x 0.71	FAD x 0.60

FAD efficiency refers to conditions at an elevation of 0 mamsl:

Temperature: 20°C
 Atmospheric pressure: 101325 Pa
 Relative humidity: 0%

6. PRODUCT DESCRIPTION

6.1. Model variations and their uses

Compressors are the source of clean, oil-free compressed air used to drive dental appliances and equipment.


Compressors models are designed for the following uses:

Dental compressors DK50 2x4VR/110 - are designed for independent installation in a suitable environment.


Dental compressors DK50 2x4VR/110/M - are designed for independent installation in a suitable environment and are equipped with an air dryer.

Dental compressors DK50 2x4VR/110S - are installed in cabinets with efficient noise mufflers and are suitable for use in office environments.


Dental compressors DK50 2x4VR/110S/M - are installed in cabinets and are equipped with air dryers.


DK50 2x4VR/110


DK50 2x4VR/110/M (NDL)


DK50 2x4VR/110S/M


DK50 2x4VR/110/M (MD)


S110R

DRYER MD WITH COOLER

DRYER NDL

COOLER

Without additional filtration equipment, the compressed air from a compressor is not suitable for the operation of breathing appliances or similar equipment.

7. FUNCTION

Compressor (Fig.1)

The compressor (1) draws in air through a filter (8) and compresses it through a check valve (3) into an air tank (2). The connected apparatus draws the compressed air from the air tank until the pressure drops to a default preset level on the air-pressure switch (4) switching the compressor on. The compressor again compresses air into the nozzle until the maximum pressure is reached and the compressor switches off. After compressor aggregate is switched off, pressure hose shall be pressure-release solenoid valve (9). Safety valve (5) prevents the pressure in air chamber from rising above the maximal allowed value. The drain valve (7) releases the condensate from the air nozzle. Compressed, dry and clean air is then ready for additional use in the air tank.

Compressor with membrane dryer (MD) (Fig. 2)

The compressor (1) draws in air through an inlet filter (8) and compresses it through a check valve (3) into a cooler (10). The air then passes through a filter (33) and micro-filter (34) into the dryer (35) with the dry and clean air passing through a check valve (3) to the air tank (2). Condensate from the filter and micro-filter is automatically released into a drainage vessel. The dryer ensures continuous drying of the compressed air. Compressed, dry and clean air is then ready for additional use in the air tank.

Compressor with NDL dryer (Fig. 3)

The compressor (1) draws in air through an inlet filter (8) and compresses it through a check valve (3) into the cooler (10). The air then proceeds through the condensate separator (15) and the drying chamber (11) with the adsorbent media where moisture is removed and then on to the air tank (2). The adsorbent media regenerates when the drying chamber is being evacuated, always after the compressor is switched off by the pressure switch. The air is then released from the adsorption chamber via the open solenoid valve as it is evacuated using dry compressed air. Air drying occurs in one chamber while regeneration occurs in the other chamber. The mode of each chamber changes in regular cycles and the drying and regeneration processes are then performed in the alternate chamber. Compressed, dry and clean air is then ready for additional use. Constant output pressure is maintained by the pressure regulator (12).

The pressure switch (13) protects the dryer from damage if the compressor operates for an extended period of time at low pressure (less than 5 bar)

Compressor cabinet (Fig. 4)

The soundproof cabinet itself is compact yet allows sufficient exchange of cooling air. The fan (18) below the compressor itself cools the compressor and operates whenever the compressor's motor is running. The cabinet's cooling fans (26) automatically start if the temperature in the cabinet reaches over 40°C if the compressor has been running for an extended period of time. The fans automatically switch off once the cabinet temperature drops to around 32°C.

An indicator (29) located on the front of the cabinet indicates maintenance is needed on the compressor with the NDL dryer unit (see the Maintenance Interval chapter)


Make sure that nothing impedes the free flow of air under and around the compressor. Never cover the hot air outlet on the top back side of the case.

If placing the compressor on a soft floor such as carpet, create space for ventilation between the base and floor or the box and floor, e.g. underpin the footings with hard pads.

Fig. 1 - Compressor

1. Compressor
2. Air tank
3. Check valve
4. Pressure switch
5. Safety valve
6. Pressure gauge
7. Condensate drain valve
8. Inlet filter
9. Solenoid valve
10. Dryer cooler
11. Adsorption dryer (NDL)
12. Pressure regulator
13. Pressure switch 2
14. Pressure relief valve
15. Condensate separator
16. Control valve
17. Check valve
18. Compressor fan
19. Box gas springs
20. Wall bump stop
21. Power supply cord
22. Pressure gauge hose
23. Angular screwing
24. Supporting bracket
25. Compressor wheel
26. Cabinet fan
27. Handle
28. Switch
29. Maintenance indicator
30. Cabinet pressure gauge
31. Temperature switch
32. Positioning truck
33. Filter
34. Micro-filter
35. Membrane dryer (MD)
36. Hour meter

Fig.2- Compressor with membrane dryer (MD)



Fig. 3 - Compressor with adsorption dryer (NDL)



Fig. 4 - Cabinet

INSTALLATION

8. USE

- The appliance must be installed and operated in a dry, well ventilated and dust-free area where ambient temperature is within the range of +5°C to +40°C and relative air humidity does not exceed 70%. Otherwise, failure-free operation of the compressor cannot be guaranteed. The compressor must be installed so that it is accessible at all times for operation and maintenance. Please ensure that the appliance label is accessible.
- The appliance must stand on a flat, sufficiently stable base. See paragraph 5 (Technical data) when positioning or lifting the compressor.
- Compressors cannot be exposed to outdoor environments. The appliance cannot be used in moist or wet environments. Do not use the compressor in the presence of explosive gases, dust or combustible liquids.
- Before connecting the compressor to medical equipment, the supplier must confirm that it meets all requirements for its use. Refer to the technical data of the product for this purpose. When a unit is to be built-in, classification and evaluation of compatibility must be done by the manufacturer or supplier of the product to be used.
- Any use other than that described in this manual is not covered by the guarantee, and the manufacturer is not liable for any damages that may result. The operator/user assumes all risk.

9. INSTALLATION

Only qualified personnel can install and start up the appliance and train operating personnel in its correct use and maintenance. Installation and training of all operators shall be confirmed by the installer's signature on the certificate of installation.

Prior to installation, ensure that the compressor is free of all transport packaging and stabilizers to avoid any risk of damage to the product.

Caution! When in operation, the compressor is hot. Burns or fire may result if contact is made by the operator or any flammable material.

Electric cord for connection to electric mains and air hoses may not be broken. The power cord may not be exposed to pulling, pressure and excessive heat.

9.1. Placement of the compressor


Fig.5 Unpacking

Dental compressor DK50 2x4VR/110 (Fig. 5)

After removing all packaging material, place the product on the floor and remove the stabilization parts (X and Y). Connect the equipment to the compressor using suitable means. Insert the mains plug into a rated outlet circuit.

Dental compressor DK50 2x4VR/110/M (MD) (Fig 5).

After unpacking the equipment, position it on the floor and free it from all packing materials and remove the fixation components (X and Y) - see Detail A. Connect the outlet pressure hose to the connector on the device. Insert the mains plug into a rated outlet circuit. Connect the hose(s) for the condensate drain to the condensate vessel.

Dental compressor DK50 2x4VR/110/M (NDL) (Fig. 5)

After removing all packaging material, place the product on the floor and remove the stabilization parts (X and Y). Place the dryer on the floor and attach it to the compressor using the bolts and washers - Fig. 7. Connect the compressor to the dryer using hoses following Fig. 10 and 11. Install the washer on the solenoid valves and connect the connectors following Fig. 9. Connect the equipment to the compressor using suitable means. Insert the mains plug into a rated outlet circuit. Connect the hose(s) for the condensate drain (from the dryer and condensate separator) to the condensate vessel.


Cabinet-mounted dental compressor DK50 2x4VR/110S (Fig.4 and Fig. 5)

After removing all packaging material, place the product on the floor and remove the stabilization parts (X and Y). Install the 2 wall bump stops on the compressor cabinet (20) on the rear and upper portion of the cabinet and position the cabinet in the desired location. The bump stops ensure that the cabinet is properly spaced from the wall for ventilation purposes. The 4 screws need to be removed and the front panel of the cabinet must be removed in order to install the compressor in the cabinet. Connect the cabinet connector (21) to the compressor. Place the pressure hose and the power cord through the opening in the rear lower portion of the cabinet following Fig. 13. Pick up the compressor using the handle and use the positioning truck (32) and the installed wheels to position it inside the cabinet. Install the cabinet's hose (22) and pressure gauge (30) into the quick coupling on the compressor and reinstall the front panel. Connect the equipment to the compressor using suitable means. Connect the mains plug to a properly rated outlet circuit.

A screwdriver must be used to disconnect the cabinet connector if the compressor is being disassembled!

(Fig.6)

Fig.6

**Cabinet-mounted dental compressor DK50 2x4VR/110S/M (MD)** (Fig.4, Fig.5)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Place the compressor into the cabinet using the procedure defined above. Before the compressor is installed in the cabinet, the hose(s) for the condensate drain must be passed through the opening on the back side of the cabinet and connected to the condensate vessel.

The condensate vessel must be on the floor in order to ensure the dryer operates properly. Otherwise the dryer may be damaged.

Cabinet-mounted compressor DK50 2x4VR/110S/M (NDL) (Fig. 4 and Fig. 5)

After unpacking the equipment, position it on the floor and remove all packing materials and the stabilization parts (X and Y) - see Detail A. Install the 2 wall bump stops on the compressor cabinet (20) on the rear and upper portion of the cabinet and position the cabinet in the desired location. The bump stops ensure that the cabinet is properly spaced from the wall for ventilation purposes. The 4 screws need to be removed and the front panel of the cabinet must be removed in order to install the compressor in the cabinet. Pick up the compressor using the handle and use the positioning truck (32) and the installed wheels to position it inside the cabinet. The pressure hoses, condensate drain hose, power cord and connectors must be pulled through the openings in the rear wall of the cabinet following Fig. 13 once the compressor is installed in the cabinet. Place the dryer on the floor and attach it to the cabinet using the bolts and washers - Fig. 8. Connect the compressor to the dryer using the pressure hoses and connectors - Fig. 13. Connect hose "1" to the dryer inlet marked "1". Connect hose "2" to the dryer outlet marked "2". Install the washer on the solenoid valves and connect the connectors following the markings - Fig. 9. Connect the condensate drain hoses from the dryer and the condensate separator to the condensate vessel - Fig. 14. Connect the cabinet connector (21) to the compressor. Install the cabinet's hose (22) and pressure gauge (30) into the quick coupling on the compressor and reinstall the front panel. Connect the equipment to the compressor using suitable means. Connect the mains plug to a properly rated outlet circuit.

The condensate vessel must be on the floor in order to ensure the dryer operates properly. Otherwise the dryer may be damaged.

Dryer installation

DK50 2x4VR/110/M

Fig.7

DK50 2x4VR/110S/M

Fig.8

Connecting the solenoid valve


Fig.9

Connecting hoses

DK50 2x4VR/110/M

Fig.10

Fig.11

Fig.12

DK50 2x4VR/110S/M

- Opening for the pressure hoses and condensate drain hose
- Opening for electrical cables and solenoid valve connectors

Fig.13

Fig.14

9.2. Compressed air outlet

(Fig.15)

Lead the pressure hose from the output of compressed air (1) to the appliance – dental set.


Fig.15

9.3. Electrical connection

Plug the electrical cord into the mains.

The appliance is equipped with a grounded plug. Make sure this connection complies with local electrical codes. The mains voltage and frequency must comply with the data stated on the appliance label.

(Fig.16)

- Keep the socket easily accessible to ensure that in an emergency the appliance can be safely disconnected from the mains.
- Connection to the power distribution box must be max.16 A.
- The connection of the earth ground pin \varnothing 6mm (1) with other appliances must be completed in accordance with local electrical codes. The female socket (2), which is not included in the standard set, is an optional accessory.

Fig.16


Electrical cable may not contact the hot parts of a compressor. Insulation could be damaged!
If any electrical cord or air hose is damaged it must be replaced immediately.

10. WIRING DIAGRAMS

3/N/PE ~ 400 V 50 Hz
 MAINS TN-S (TN-C-S)
 ELECTRIC OBJECT 1ST.CAT
 TYPE B


E1 - E4	Fan of compressor	M1, M2	Motor of compressor
E5 - E6	Fan of cooler	F1, F2	Breaker
M12	Solenoid valve	K3	PCB
X1	Terminal with fuse	Q1, Q2	Contact
F3	Fuse	B2	Thermo switch
X2	Connector	B1	Pressure switch
		P1	Hour meter

DK50 2x4VR/110/M (MD)

DK50 2x4VR/110


3/N/PE ~ 400/230 V 50.60 Hz
 MAINS TN-S [TN-C-S]
 ELECTRIC OBJECT 1ST.CAT
 TYPE B

DK50 2x4VR/110


- | | | | |
|---------|-------------------------|---------|---------------------|
| E1 - E4 | Fan of compressor | M1,M2 | Motor of compressor |
| E5 - E8 | Fan of cooler | F1,F2 | Breaker |
| M12,M13 | Solenoid valve | K3 | PCB |
| M10,M11 | Solenoid valve of dryer | Q11,Q12 | Contactur |
| K1 | Control unit | B2 | Thermo switch |
| X1 | Terminal with fuse | B1,B3 | Pressure switch |
| B4 | Reed switch | P1 | Hour meter |
| F3 | Fuse | | |
| X2 | Connector | | |

3/N/PE ~ 400 V 50 Hz
 MAINS TN-S [TN-C-S]
 ELECTRIC OBJECT 1ST.CAT
 TYPE B

DK50 2x4VR/110/M(NDL)


3/N/PE - 400 V 50 Hz
 MAINS TN-S [TN-C-S]
 ELECTRIC OBJECT 1ST.CAT
 TYPE B

DK50 2x4VR/110S/M(NDL)

1/N/PE - 230V 50/60Hz
 ELEKTRIC OBJECT 1st.CAT
 TYPE B

S110-4VR


- | | |
|----------|----------------------------|
| M1, M2 | Motor of compressor |
| Q11, Q12 | Contactors |
| FA1, FA2 | Breaker |
| B2 | Thermo switch |
| B1, B3 | Pressure switch |
| E1 - E4 | Fan of compressor |
| E5 - E8 | Fan of cooler |
| M12 | Solenoid valve |
| M10, M11 | Solenoid valve of dryer |
| K1 | Control unit |
| X1 | Terminal with fuse |
| B4 | Reed switch |
| F3 | Fuse |
| E10-E13 | Fan of box |
| X2 | Connector |
| X10 | Plug |
| S10 | Switch |
| P10 | LED-SERVIS for NDL |
| B10 | Cabinet temperature switch |
| P1 | Hour meter |


11. PNEUMATIC DIAGRAM

- 1 Inlet filter
- 2 Compressor
- 3 Fan
- 4 Solenoid valve
- 5 Safety valve
- 6 Check valve
- 7 Pressure switch
- 8 Air tank
- 9 Pressure gauge
- 10 Dryer outlet solenoid valve
- 11 Cooler
- 12 Cooler fan
- 13 Condensate separator
- 14 Pressure relief valve
- 15 Pressure regulator
- 16 Dryer
- 17 Vessel
- 18 Filter
- 19 Mikrofilter


DK50 2x4VR/110

DK50 2x4VR/110/M (NDL)

DK50 2x4VR/110/M (MD)

12. FIRST OPERATION

(Fig.17)

- Make sure that all stabilizers used during transport were removed.
- Check that all pressurized air line connections are secure.
- Connect to the mains.
- Start compressor at pressure switch (2) by turning switch (3) to position "I."
- For cabinet-mounted compressors, turn the switch (28) in Fig. 4 on the front side of the equipment cabinet into the "I" position; a lighted switch indicates the equipment is in operation.

• **Compressor** - At first operation the air tank is pressurized until it reaches a preset level when the compressor automatically switches off. As the air is used, the compressor works in automatic mode, switched on or off by the pressure switch.

• **Compressor with dryer** – Compressor operation is the same and the moisture from the compressed air is removed as it passes through the dryer.

For NDL - the condensate is also blown out of the drain on the dryer, which can be heard as a brief whoosh sound when the compressor shuts down or when the drying chambers switch.

The compressor is not equipped with an emergency power supply.

OPERATION

In case of emergency, disconnect the compressor from the mains (pull out the mains plug).

**The compressor has hot surfaces.
Burns or fire may result if contact is made.**

During prolonged operation of the compressor, the temperature in the box may increase to over 40°C. At this point the cooling fan automatically switches on. After cooling the space to under 32°C, the ventilator switches off.

Automatic start: when pressure in the tank drops to the pressure switch's lower limit level, the compressor automatically switches on. The compressor automatically switches off after reaching the pressure switch's upper limit level.

Compressor with dryer

A correct function of the drier depends on the compressor's operation and no attendance is required. The pressure vessel need not be sludged, because the pressure air entering the air chamber is already dried.

- It is forbidden to alter the working pressures of pressure switch set by manufacturer. The operation of the compressor at working pressure lower than the switching pressure demonstrates the overload of the compressor (high air consumption) by the appliance, leakages in pneumatic distributions, failure of aggregate or drier.
- Prior connecting drier to air chamber, that was used with compressor without drier, it is necessary to clean interior surface of air chamber and perfectly remove condensed liquid. Then interconnect electric part of drier with compressor according to wiring diagram in accord with valid regional regulations.

Required drying performance can only be achieved when following the defined operating conditions!

Drying performance will decline and the achieved dew point will drop if the dryer is operated at any pressure below the minimum working pressure!
Dryer operation at a pressure of 0.5 Bar below the minimum working pressure can lower the dew point at the outlet by more than 10°C!

The dryer will be irrevocably damaged and need replacement if operated at any temperature above the maximum working temperature

13. SWITCHING THE COMPRESSOR ON

(Fig.17)

Switch on the compressor at the pressure switch (2) by turning the knob (3) to position "I." (for compressor in the box switch (28) Fig.4), on the front part of the compressor box), The compressor sends pressurized air to the air tank. As the compressed air is used, the pressure in the air nozzle drops to a preset level, the compressor switches on and the air nozzle fills with compressed air. After reaching the cutoff pressure the compressor turns off automatically and the cycle is repeated. Check the value of switching-on and switching-off pressure on pressure gauge. The values may be within a tolerance of $\pm 10\%$. Air pressure in air chamber must not exceed maximal permitted operation pressure..



Fig.17

Never tamper with the pressure switch (2). Adjustments are not allowed. The pressure switch (2) has been set by the manufacturer and further setting of switching on and off pressure may be carried out only by a qualified expert trained by the manufacturer.


CONTROL PANEL SCREENS FOR THE NDL DRYER

- Supplier
- Supplier contact
- Operating hours
- Chamber A (B) in operation


Status messages

- Low pressure


- a) From compressor start-up until pressure reaches 5.5 bar
- b) During operation if pressure drops below 5.1 bar

- 24-hour regeneration mode


Dryer regeneration for 10 minutes is automatically switched on after 24 hours of compressor operation.

- 2 year maintenance

Calculated when equipment is switched on

- 8,000 hour maintenance

Calculated using the number of operating hours

MAINTENANCE

14. MAINTENANCE SCHEDULE

Notice!

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used.

Time interval	Maintenance that must be performed	Chapter	Performed by
1 x day	Emptying condensate At high air humidity	15.2	User
1x per week 1x per week - check function	Compressors without air dryers Compressors with air dryers		
Once annually	Equipment cleaning	15.1	Qualified technician
	Check of dryer function	15.3	
	Safety valve check	15.4	
	NDL Replace float in water separator	15.6	
	MD Replace filter and micro-filter elements	15.7 15.8	
	Check tightness of joints and overall inspection of the equipment	Maintenance documentation	
1 x every 2 years	Conduct "repeated test" pursuant to EN 62353	14	Qualified technician
1 x every 2 years or after 5000 hours	Replace inlet filter and pre-filter	15.5	
8,000 hours (or 2 years) 16,000 hours (or 4 years) 24,000 hours (or 6 years)	NDL Replacement maintenance kit - dryer	15.9	Qualified technician
	Replace media in both dryer cylinders and seals.		
	Replace media in both dryer cylinders and seals.		
	Replace media in both dryer cylinders, seals and all valves.		

15. MAINTENANCE

Repair work beyond normal maintenance can be performed only by qualified personnel or the manufacturer's representative.
Use only spareparts and accessories approved by the manufacturer.

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).

TO ENSURE THAT THE COMPRESSOR WORKS CORRECTLY, PERFORM THE FOLLOWING MAINTENANCE TASKS AT REGULAR INTERVALS (CHAPTER 14):

THE CABINET MUST BE OPENED ON CABINET-MOUNTED COMPRESSORS BEFORE CONDUCTING THE FOLLOWING CHECKS.

15.1. Equipment cleaning

The equipment, in particular the cooling fan must be kept clean to ensure the long-term efficiency of the cooler – the cooling fins must be cleaned occasionally to remove any dust.

15.2. Condensation drain valve

Compressors (Fig.18)

Switch off the compressor at the mains. Reduce air pressure in the appliance to max. 1 bar by releasing air via a connected device. Place the vessel under release valve (1) and open it. Wait until condensation is fully drained from the pressure tank. Close drain valve (1).


Fig.18

Compressors with air dryer

In the case of a regular operation condensate is automatically excreted via air dryer and it is entrapped in a bottle. Take out the bottle from a holder, release cap and pour out the condensate.

15.3. Check of dryer function

Check for proper dryer function by opening the vent valve. No condensate should leak out. Seek professional assistance if condensate appears during such venting!

15.4. Safety valve check

(Fig.17)

When the compressor is operated for the first time, make sure that the safety valve is working properly. Turn screw (4) of safety valve (1) several rotations to the left until the safety valve releases air. Let the safety valve blow out for only a few seconds. Turn screw (4) to the right until it seats, closing the valve.

The safety valve must never be used for depressurizing the air tank. It could damage the safety valve. The valve is set to the maximum permitted pressure by the manufacturer. Adjustments are not permitted.

Warning! Compressed air can be dangerous. Wear eye protection when blowing air out.

15.5. Replacement of the inlet filter and prefilter

(Fig.19)


At the lid of the compressors crankcase is an Inlet filter (1) and prefilter (3).

Replacing of the inlet filter:

- Hand pull the rubber stopper (2).
- Remove used and dirty filter.
- Input new filter and set rubber stopper.

Replacing of the prefilter:

- Hand pull prefilter (3).
- Replace old prefilter with new

	Order number
Input filter	025200139
Prefilter	025200150

15.6. Replace float in water separator

(Fig.20)

During regular operation of drier is necessary to change the float in the water separator.

- Release pressure.
- Dismount separator container.
- Pull out condenser separator
- Release nut of the float placed in the bottom part of the container.
- Pull out worn float of the separator and replace it for the new one.
- Secure the float with the nut in bottom part of the container.
- Place the condenser separator back as shown in the picture.
- Container of the separator insert back and screw in.
- The container is definitely secured in the point of the symbol.


Fig.20

Water separator	Order number	Float	Order number
WS010BBFX	025200119	EF1	025200146

15.7. Replace the filter element

(Fig. 21)

Release the catch (1) on the filter vessel by pulling down.

Rotate the vessel (2) and remove.

Unscrew the filter retainer (3).

Replace the filter element (4) and reinstall the filter retainer.

Reinstall the filter vessel and rotate to secure against the catch.

Fig. 21

Filter	Order number	Filter element	Order number
AF 40P-060S	025200215-000	AF 40P-060S 5 μm	025200079-000

15.8. Replacing the micro-filter element

(Fig. 22)

Release the catch (1) on the micro-filter vessel by pulling down.

Rotate the vessel (2) and remove.

Remove the filter (3).

Install the new filter element.

Reinstall the filter vessel and rotate to secure against the catch.


Fig. 22

Micro-filter	Order number	Filter element	Order number
AFM 40-F03C	0252002156-000	AFM 40P-060AS 5 μm	025200080-000

15.9. NDL dryer maintenance kit replacement

Before starting any work, disconnect the equipment from all power sources and vent the air tank to zero

Repair work beyond normal maintenance can only be performed by qualified personnel or the manufacturer's representative.

The maintenance kit contains: Cartridge with drying media, gaskets and valves depending on maintenance interval.

An indicator indicates maintenance interval :

- on the display - Zelio Logic
- on the cabinet


Fig.23-A: Spare parts


Fig.23-B: Dryer cartridge

Dryer	Order number	Dryer service kit	Order number
NDL-60	035900116	NDK-30	

16. STORAGE

If the compressor will not be used for a prolonged time period, drain any condensate from the air tank. Then turn on the compressor for 10 minutes, keeping the drain valve open (1) (Fig.18). Switch off the compressor by switch (3) at pressure switch (2) (Fig.17), close the drain valve and disconnect the appliance from the mains.

17. DISPOSING OF THE APPLIANCE

- Disconnect the appliance from the mains.
- Release air pressure in the pressure tank by opening the drain valve (1) (Fig.18).
- The components of the product are non-toxic.
- Dispose of the appliance following all environmental regulations.

18. REPAIR SERVICE

Guaranteed and post-guarantee repairs must be done by the manufacturer, its authorized representative, or service personnel approved by the supplier.

The manufacturer reserves the right to make changes to the appliance without notice. Any changes made will not affect the functional properties of the appliance.

19. SOLVING PROBLEMS

Caution! Before proceeding, depressurize the air tank to zero and disconnect the appliance from the mains.

The internal surfaces of the air tank must be cleaned and all condensed liquid must be removed after a dryer failure.

Check the dew point of the air leaving the air tank (see Chapter 5 - Technical Data) in order to protect connected equipment from damage!

Troubleshooting can be performed only by qualified personnel.

FAILURE	POSSIBLE CAUSE	REMEDY
Compressor does not start	No voltage in pressure switch Disconnected winding of motor, damaged thermal protection Faulty capacitor Seizure of piston or another rotary part Pressure switch does not switch on	Check voltage in socket Check fuse – replace faulty one Loosen terminal – tighten it Check power cord – replace faulty one Replace motor or re-wind it Replace capacitor Replace damaged parts Check the function of pressure switch
Compressor often switches on	Air leak in pneumatic distribution system Leaking check valve Greater volume of condensed liquid in pressure vessel	Check pneumatic distribution system – seal loose joint Clean valve, replace seals, replace valve Drain condensed liquid
Prolonged running of compressor	Air leak in pneumatic distribution system Worn piston ring Contaminated input filter and prefilter Dirty filter in the dryer Defective solenoid valve	Check pneumatic distribution system – seal loose joint Replace worn piston ring Replace contaminated filters with the new ones Change the outlet filter and inspect dessicant Repair or change the valve
Compressor is noisy (knocking, metal noises)	Damaged bearing of piston, piston rod, motor bearing Loose or cracked spring	Replace damaged bearing Replace damaged spring
	Membrane dryer	
Dryer doesn't dry (condensed water in the tank)	inoperative cooler ventilator	replace ventilator ,check supply of electric energy
	Damaged dryer	Replace dryer
	Dirty automatic condensate drain on filters Dirty filter and micro-filter elements	clean / replace Replace old elements with new elements
	NDL Dryer	
Decreased dew point performance	1. Insufficient pressure at inlet 2. Electrical failure 3. Damp or contaminated dryer media 4. Excessively high air consumption 5. Excessive inlet air temperature 6. Insufficient clean air 7. Exhaust silencer blocked	1. Inlet pressure must be at least 4 Bar. If this is not the case, adjust and set the inlet pressure. 2. Ensure the equipment is switched on and the dryer's front panel is lit, check for proper dryer cycle. 3. Remove the source of the contamination. Replace the cartridges – do not reuse. 4. Ensure dryer performance matches required air demand. 5. Check technical specifications. 6. Cancel incorrect settings and contact a service technician – adjust settings. 7. Contact a service technician.
Dryer cycle failure	8. Controller malfunction 9. Lamp is not lit 10. Insufficient inlet pressure 11. Venting malfunction during regeneration 12. Output flow halted	8. Check to ensure controller has power, check the display to ensure the solenoid valves are powered in regular cyclical operation. 9. Check unit power and fuses 10. Inlet pressure must be at least 4 Bar. If this is not the case, adjust and set the inlet pressure. 11. If the solenoid valve has power and is malfunctioning – replace the valve. The valve is working properly if a click can be heard at the outlet when venting 12. Check the air flow at the inlet.
Constant venting	13. Dryer run failure 14. Irregular air flow from exhaust	13. Switch off and restart the dryer. Check to see if the dryer is under pressure before switching on – to ensure the dryer can be started before starting operation. 14. Defective or damaged valve, maintenance required

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