



Model: DA-100

- * Image Sensor: CMOS ¼
- * Camera Resolution: 15 megapixels
- * OS: Embedded Operating System
- * Screen Size: 17 inch
- * Native Resolution: 1280*1024
- * Picture Display: 1 Picture/4 Pictures
- * Storage: USB Flash Device (Max 32GB)
- * Picture Storage Format: JPG
- * Focus Range: 5mm-50mm
- * WIFI Transmission:
 - Scan QR code to get images
 - Real-time operation with PC
- * Other Functions:
 - Video Recording
 - Multimedia
- * Color Option: White, Beige

Keys Functions:

- ① Image Capture, Save Images, Turn off the camera.
- ② Switch Images Display, Video Recording.
- ③ Page Up, Create QR Code.
- ④ Page Down
- ⑤ Cancel the image, Delete single image, Delete all the images.

